


SCGN Antibody / Secretagogin [clone SCGN/7326] (V5160)


Catalog No.	Formulation	Size
V5160-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V5160-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V5160SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	SCGN/7326
Purity	Protein A/G affinity
UniProt	O76038
Localization	Cytoplasm, Secreted
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This SCGN antibody is available for research use only.



IHC staining of FFPE human brain tissue with SCGN antibody (clone SCGN/7326). Inset: PBS used in place of primary Ab (secondary Ab negative control). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free SCGN antibody (clone SCGN/7326) as confirmation of integrity and purity.

Description

SCGN, also known as Secretagoin, CALBL, setagin or SECRET, is a 276 amino acid cytoplasmic protein that contains six EF-hand domains and is related to the calcium-binding proteins Calretinin and Calbindin D28K. Expressed in a variety of tissues including stomach, thyroid, colon, brain and neuroendocrine cells, SCGN is thought to be involved in cell proliferation and KCl (potassium chloride)-mediated calcium flux events. Through its interaction with KCl and its subsequent ability to modulate calcium storage pools within the cell, SCGN may function to negatively control growth and differentiation rates and, thus, indirectly inhibit cell replication.

Application Notes

Optimal dilution of the SCGN antibody should be determined by the researcher.

Immunogen

A recombinant partial protein sequence (within amino acids 100-276) from the human protein was used as the immunogen for the SCGN antibody.

Storage

Aliquot the SCGN antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.