


SCOC Antibody (R34305)

Catalog No.	Formulation	Size
R34305-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	60592
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SCOC antibody is available for research use only.



Description

Additional name(s) for this target protein: Short coiled-coil protein

Application Notes

Optimal dilution of the SCOC antibody should be determined by the researcher.

Immunogen

Amino acids HTLEDLSARVDAVKE were used as the immunogen for this SCOC antibody. This sequence is common to isoforms 1-6. **Storage** Aliquot and store the SCOC antibody at -20oC.