


Scraps Antibody (R36348)

Catalog No.	Formulation	Size
R36348-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54443
Applications	Immunohistochemistry (FFPE) : 3-10ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Scraps antibody is available for research use only.



IHC testing of FFPE human kidney with Scraps antibody at 10ug/ml. HIER: Microwave with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: Anillin; ANLN

Application Notes

Optimal dilution of the Scraps antibody should be determined by the researcher.

Immunogen

Amino acids WQPDACYKPIGKP were used as the immunogen for this Scraps antibody.

Storage

Aliquot and store the Scraps antibody at -20oC.