


SDCCAG8 Antibody (R36187)

Catalog No.	Formulation	Size
R36187-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10806
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SDCCAG8 antibody is available for research use only.



Western blot testing of HEK293 lysate with SDCCAG8 antibody at 0.3ug/ml. Predicted molecular weight: ~82kDa.

Description

Additional name(s) for this target protein: Serologically defined colon cancer antigen 8

Application Notes

Optimal dilution of the SDCCAG8 antibody should be determined by the researcher.

Immunogen

Amino acids TINDQSQYIHHLEAE were used as the immunogen for this SDCCAG8 antibody.

Storage

Aliquot and store the SDCCAG8 antibody at -20oC.