

SECISBP2 Antibody (R34416)


Catalog No.	Formulation	Size
R34416-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79048
Applications	Western Blot : 2-6ug/ml Immunohistochemistry (FFPE) : 3ug/ml ICC/IF : suitable Immunoprecipitation : suitable ELISA (peptide) LOD : 1:64000
Limitations	This SECISBP2 antibody is available for research use only.



Western blot testing of HFF cell lysate and A) control siRNA or B) 100nM siRNA with SECISBP2 antibody at 0.2ug/ml.



IHC staining of FFPE human testis with SECISBP2 antibody at 3ug/ml. HIER: microwaved with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: SECIS binding protein 2, SBP2

Application Notes

Optimal dilution of the SECISBP2 antibody should be determined by the researcher.

Immunogen

Amino acids KERQERKQRLQEN were used as the immunogen for this SECISBP2 antibody.

Storage

Aliquot and store the SECISBP2 antibody at -20oC.