


Secretogranin II Antibody (R35245)

Catalog No.	Formulation	Size
R35245-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7857
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Secretogranin II antibody is available for research use only.



Western blot testing of human cerebellum lysate with Secretogranin II antibody at 0.3ug/ml. Predicted molecular weight: ~71kDa.

Description

Additional name(s) for this target protein: SCG2, CHGC, Chromogranin-C

Application Notes

Optimal dilution of the Secretogranin II antibody should be determined by the researcher.

Immunogen

Amino acids QELGKLTGPNNQKRE were used as the immunogen for this Secretogranin II antibody.

Storage

Aliquot and store the Secretogranin II antibody at -20oC.