


# SENP6 Antibody (R34569)

Catalog No.	Formulation	Size
R34569-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26054
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SENP6 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: SUMO1/sentrin specific peptidase 6

## **Application Notes**

Optimal dilution of the SENP6 antibody should be determined by the researcher.

### **Immunogen**

Amino acids KPKYEPNPHYHEN were used as the immunogen for this SENP6 antibody. **Storage** Aliquot and store the SENP6 antibody at -20oC.