

SENP6 Antibody (R34569)

Catalog No.	Formulation	Size
R34569-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26054
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SENP6 antibody is available for research use only.



Western blot testing of HeLa cell nuclear lysate with SENP6 antibody at 0.5ug/ml.
Predicted molecular weight: ~125kDa, observed here at ~150kDa.

Description

Additional name(s) for this target protein: SUMO1/sentrin specific peptidase 6

Application Notes

Optimal dilution of the SENP6 antibody should be determined by the researcher.

Immunogen

Amino acids KPKYEPNPHYHEN were used as the immunogen for this SENP6 antibody.

Storage

Aliquot and store the SENP6 antibody at -20oC.