


SEPT4 Antibody (R33875)

Catalog No.	Formulation	Size
R33875-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5414
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This SEPT4 antibody is available for research use only.



Western blot testing of human cerebellum lysate with SEPT4 antibody at 0.5ug/ml.
Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: Septin 4

Application Notes

Optimal dilution of the SEPT4 antibody should be determined by the researcher.

Immunogen

Amino acids ELSKFVKDFSGN were used as the immunogen for this SEPT4 antibody.

Storage

Aliquot and store the SEPT4 antibody at -20°C.