


Septin 4 Antibody (R35692)

Catalog No.	Formulation	Size
R35692-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5414
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Septin 4 antibody is available for research use only.



Western blot testing of human heart lysate with Septin 4 antibody at 1ug/ml. Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: ARTS (Apoptosis-related protein in the TGF-beta signaling pathway), BRADEION, CE5B3, H5; SEPT4

Application Notes

Optimal dilution of the Septin 4 antibody should be determined by the researcher.

1. This Septin 4 antibody will detect isoform 1/3/4.

Immunogen

Amino acids RNKLTRESGTD were used as the immunogen for this Septin 4 antibody.

Storage

Aliquot and store the Septin 4 antibody at -20oC.