


Septin 4 Antibody (R35692)

Catalog No.	Formulation	Size
R35692-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5414
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Septin 4 antibody is available for research use only.



Description

Additional name(s) for this target protein: ARTS (Apoptosis-related protein in the TGF-beta signaling pathway), BRADEION, CE5B3, H5; SEPT4

Application Notes

Optimal dilution of the Septin 4 antibody should be determined by the researcher.

1. This Septin 4 antibody will detect isoform 1/3/4.

Immunogen

Amino acids RNKLTRESGTD were used as the immunogen for this Septin 4 antibody.

Storage

Aliquot and store the Septin 4 antibody at -20oC.