


Septin 7 Antibody (R33560)

Catalog No.	Formulation	Size
R33560-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Horse, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	989
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Septin 7 antibody is available for research use only.



Description

Additional name(s) for this target protein: SEPT7; CDC10; CDC3

Application Notes

Optimal dilution of the Septin 7 antibody should be determined by the researcher.

Immunogen

Amino acids DNNKNKGQLTKSP were used as the immunogen for this Septin 7 antibody. **Storage** Aliquot and store the Septin 7 antibody at -20oC.