


Septin 7 Antibody (R33560)

Catalog No.	Formulation	Size
R33560-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Horse, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	989
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Septin 7 antibody is available for research use only.



Western blot testing of nuclear HeLa lysate with Septin 7 antibody at 0.5ug/ml. Predicted molecular weight: ~51kDa.

Description

Additional name(s) for this target protein: SEPT7; CDC10; CDC3

Application Notes

Optimal dilution of the Septin 7 antibody should be determined by the researcher.

Immunogen

Amino acids DNNKNKGQLTKSP were used as the immunogen for this Septin 7 antibody.

Storage

Aliquot and store the Septin 7 antibody at -20oC.