


SET Antibody (isoform 2) (R35837)


Catalog No.	Formulation	Size
R35837-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6418
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SET antibody is available for research use only.



Western blot testing of human kidney lysate with SET antibody at 0.3ug/ml. Predicted molecular weight: 32kDa.



Western blot testing of Daudi (A) and MOLT4 (B) cell lysate using SET antibody at 0.01ug/ml. Predicted molecular weight: 32kDa.

Description

Additional name(s) for this target protein: SET nuclear oncogene, I2 alpha PP2A, 2PP2A, I2PP2A, IGAAD, IPP2A2

Application Notes

Optimal dilution of the SET antibody should be determined by the researcher.

1. This SET antibody is specific to isoform 2.

Immunogen

Amino acids SAPAAKVSKKELNS were used as the immunogen for this SET antibody. This sequence is specific to isoform 2.

Storage

Aliquot and store the SET antibody at -20oC.