


SETMAR Antibody (R35812)

Catalog No.	Formulation	Size
R35812-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6419
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This SETMAR antibody is available for research use only.



Western blot testing of Daudi lysate with SETMAR antibody at 2ug/ml. Predicted molecular weight: ~78kDa.

Description

Additional name(s) for this target protein: SET domain and mariner transposase fusion

Application Notes

Optimal dilution of the SETMAR antibody should be determined by the researcher.

Immunogen

Amino acids RWQKCVDCNGSYFD were used as the immunogen for this SETMAR antibody.

Storage

Aliquot and store the SETMAR antibody at -20oC.