


SGSM3 Antibody (R35558)

Catalog No.	Formulation	Size
R35558-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27352
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This SGSM3 antibody is available for research use only.



Western blot testing of A549 lysate with SGSM3 antibody at 2ug/ml. Predicted/observed molecular weight: ~85kDa.

Description

Additional name(s) for this target protein: Small G protein signaling modulator 3, RUTBC3, RUN and SH3 containing 3

Application Notes

Optimal dilution of the SGSM3 antibody should be determined by the researcher.

Immunogen

Amino acids RDHENYVACSRSHRR were used as the immunogen for this SGSM3 antibody.

Storage

Aliquot and store the SGSM3 antibody at -20oC.