


SH3BP1 Antibody (R35121)


Catalog No.	Formulation	Size
R35121-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23616
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 0.3-1ug/ml Immunocytochemistry : suitable ELISA (peptide) LOD : 1:16000
Limitations	This SH3BP1 antibody is available for research use only.



Western blot testing of human brain lysate with SH3BP1 antibody at 1ug/ml. Predicted molecular weight: ~76kDa.



IHC staining of FFPE human kidney with SH3BP1 antibody at 0.3ug/ml. Figure A shows cytoplasmic staining of lymphocytes in the interstitial parenchyma of the kidney. Figure B shows staining in a glomerulus of the kidney. HIER: microwaved with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: SH3-domain binding protein 1

Application Notes

Optimal dilution of the SH3BP1 antibody should be determined by the researcher.

Immunogen

Amino acids PQPRPRSLASETN were used as the immunogen for this SH3BP1 antibody.

Storage

Aliquot and store the SH3BP1 antibody at -20oC.