

## SHP-1 Antibody / PTPN6 (R34303)

Catalog No.	Formulation	Size
R34303-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	5777
<b>Applications</b>	Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This SHP-1 antibody is available for research use only.



Western blot testing of Jurkat cell lysate with SHP-1 antibody at 0.2ug/ml. Predicted molecular weight: ~68kDa.

### Description

Additional name(s) for this target protein: Protein tyrosine phosphatase, non-receptor type 6, HCP, Hematopoietic cell phosphatase; PTPN6

### Application Notes

Optimal dilution of the SHP-1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids HTKNKREEKVKKQ were used as the immunogen for this SHP-1 antibody.

## **Storage**

Aliquot and store the SHP-1 antibody at -20oC.