


SIAH1 Antibody Biotin Conjugate (R36054BTN)

Catalog No.	Formulation	Size
R36054BTN	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6477
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This SIAH1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Seven in absentia homolog 1

Application Notes

Optimal dilution of the SIAH1 antibody should be determined by the researcher.

Immunogen

Amino acids SRQTATALPTGTSKC were used as the immunogen for this SIAH1 antibody. **Storage** Aliquot and store the SIAH1 antibody at -20oC.