

SIAH1 Antibody Biotin Conjugate (R36054BTN)

| Catalog No. | Formulation | Size |
|-------------|---|--------|
| R36054BTN | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse, Rat |
| Predicted Reactivity | Cow, Dog |
| Format | Biotin Conjugate |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 6477 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This SIAH1 antibody is available for research use only. |



Western blot testing of mouse brain lysate with biotinylated SIAH1 antibody at 1 ug/ml.
Predicted molecular weight: ~31kDa.

Description

Additional name(s) for this target protein: Seven in absentia homolog 1

Application Notes

Optimal dilution of the SIAH1 antibody should be determined by the researcher.

Immunogen

Amino acids SRQTATALPTGTSKC were used as the immunogen for this SIAH1 antibody.

Storage

Aliquot and store the SIAH1 antibody at -20oC.