


SIAH1 Antibody / Seven in absentia homolog 1 (R36054)


Catalog No.	Formulation	Size
R36054-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6477
Applications	Western Blot : 1-3ug/ml Immunoprecipitation : suitable ELISA (peptide) LOD : 1:16000
Limitations	This SIAH1 antibody is available for research use only.



Western blot testing of human liver lysate with SIAH1 antibody at 1ug/ml. Predicted molecular weight: ~31kDa.



Western blot of A) mouse brain, B) rat brain, C) mouse liver and D) rat liver lysates with SIAH1 antibody at 1ug/ml. Predicted molecular weight: ~31kDa.

Description

Additional name(s) for this target protein: Seven in absentia homolog 1

Application Notes

Optimal dilution of the SIAH1 antibody should be determined by the researcher.

Immunogen

Amino acids SRQTATALPTGTSKC were used as the immunogen for this SIAH1 antibody.

Storage

Aliquot and store the SIAH1 antibody at -20oC.