


SIPA1 Antibody (R33174)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33174-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|---|
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 6494 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This SIPA1 antibody is available for research use only. |



Western blot testing of human EBV-immortilized lymphoblastoid lysate with SIPA1 antibody at 1ug/ml, with (B) and without (A) immunizing peptide. Predicted molecular weight: ~112kDa but routinely observed at ~130kDa. It is observed here at ~150kDa.

Description

Additional name(s) for this target protein: Signal-induced proliferation-associated gene 1

Application Notes

Optimal dilution of the SIPA1 antibody should be determined by the researcher.

Immunogen

Amino acids ASKQLGSPTADLA were used as the immunogen for this SIPA1 antibody.

Storage

Aliquot and store the SIPA1 antibody at -20oC.

References (1)