


SIPA1 Antibody (R33174)

Catalog No.	Formulation	Size
R33174-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6494
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SIPA1 antibody is available for research use only.



Western blot testing of human EBV-immortalized lymphoblastoid lysate with SIPA1 antibody at 1ug/ml, with (B) and without (A) immunizing peptide. Predicted molecular weight: ~112kDa but routinely observed at ~130kDa. It is observed here at ~150kDa.

Description

Additional name(s) for this target protein: Signal-induced proliferation-associated gene 1

Application Notes

Optimal dilution of the SIPA1 antibody should be determined by the researcher.

Immunogen

Amino acids ASKQLGSPTADLA were used as the immunogen for this SIPA1 antibody.

Storage

Aliquot and store the SIPA1 antibody at -20oC.

References (1)