


## SLAMF8 Antibody (R34154)

Catalog No.	Formulation	Size
R34154-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	56833
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SLAMF8 antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human SLAMF8-MYC probed with SLAMF8 antibody [0.01ug/ml] in (A) and anti-MYC [1/1000] in (C). Mock-transfected HEK293 probed with SLAMF8 antibody [0.01ug/ml] in (B). Predicted molecular weight: ~32/45kDa (unmodified/glycosylated).

## Description

Additional name(s) for this target protein: SLAM family member 8; BLAME; CD353

## Application Notes

Optimal dilution of the SLAMF8 antibody should be determined by the researcher.

## Immunogen

Amino acids HADRVGPETENP were used as the immunogen for this SLAMF8 antibody.

## Storage

Aliquot and store the SLAMF8 antibody at -20oC.