


SLC17A5 Antibody (R33800)

Catalog No.	Formulation	Size
R33800-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9NRA2
Gene ID	26503
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SLC17A5 antibody is available for research use only.



Western blot testing of human placenta lysate with SLC17A5 antibody at 1ug/ml. Predicted molecular weight ~55 kDa and 31 kDa (two isoform) which may be observed at higher molecular weights due to glycosylation.

Description

Additional name(s) for this target protein: Solute carrier family 17; Sialin; HP59

Application Notes

Optimal dilution of the SLC17A5 antibody should be determined by the researcher.

Immunogen

Amino acids EHSAPIKVHHNQTGK were used as the immunogen for this SLC17A5 antibody. This sequence is common to isoforms 1 and 2.

Storage

Aliquot and store the SLC17A5 antibody at -20oC.