


## SLC22A3 Antibody / OCT3 (R35997)

Catalog No.	Formulation	Size
R35997-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6581
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This SLC22A3 antibody is available for research use only.



Western blot testing of HEK293 lysate with SLC22A3 antibody at 0.5ug/ml. Predicted molecular weight: ~61kDa.

### Description

Additional name(s) for this target protein: Solute carrier family 22 (extraneuronal monoamine transporter), member 3; OCT3 (organic cation transporter 3)

### Application Notes

Optimal dilution of the SLC22A3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SNYSEITVTDEEVSN were used as the immunogen for this SLC22A3 antibody.

## **Storage**

Aliquot and store the SLC22A3 antibody at -20oC.