


SLC35A4 Antibody / Solute carrier family 35 member A4 (RQ7760)

Catalog No.	Formulation	Size
RQ7760	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q96G79
Applications	Western Blot : 0.5-1ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This SLC35A4 antibody is available for research use only.



Western blot testing of 1) human Caco-2, 2) human A549, 3) rat thymus and 4) mouse thymus tissue lysate with SLC35A4 antibody. Predicted molecular weight ~35 kDa.

Description

Predicted to enable pyrimidine nucleotide-sugar transmembrane transporter activity. Involved in positive regulation of translation in response to stress. Predicted to be located in Golgi apparatus. Predicted to be integral component of Golgi membrane.

Application Notes

Optimal dilution of the SLC35A4 antibody should be determined by the researcher.

Immunogen

E. coli-derived recombinant human protein (amino acids M1-Q142) was used as the immunogen for the SLC35A4 antibody.

Storage

After reconstitution, the SLC35A4 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.