


# SLC40A1 Antibody (R33865)

Catalog No.	Formulation	Size
R33865-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	30061
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SLC40A1 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Solute carrier family 40 (iron-regulated transporter), member 1; Ferroportin 1; FPN; SLC40A1

#### **Application Notes**

Optimal dilution of the SLC40A1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids ELKQLNLHKDTEPKP were used as the immunogen for this SLC40A1 antibody.

# **Storage**

Aliquot and store the SLC40A1 antibody at -20oC.