


## SLC47A2 Antibody (R35865)

Catalog No.	Formulation	Size
R35865-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	146802
Applications	Western Blot : 0.1-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SLC47A2 antibody is available for research use only.



Western blot testing of human testis lysate with SLC47A2 antibody at 0.3ug/ml.  
Predicted molecular weight: ~60kDa.

## Description

Additional name(s) for this target protein: Solute carrier family 47, member 2

## Application Notes

Optimal dilution of the SLC47A2 antibody should be determined by the researcher.

## Immunogen

Amino acids SECHVDFFRTPEEAH were used as the immunogen for this SLC47A2 antibody.

## **Storage**

Aliquot and store the SLC47A2 antibody at -20oC.