


## SNAI1 Antibody (F47916)


Catalog No.	Formulation	Size
F47916-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47916-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	O95863
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This SNAI1 antibody is available for research use only.



Western blot analysis of lysate from (1) HeLa and (2) SiHa cell line using SNAI1 antibody at 1:1000. Predicted molecular weight ~29kDa.



Western blot analysis of mouse kidney lysate using SNAI1 antibody at 1:1000.

## Description

The Drosophila embryonic protein snail is a zinc finger transcriptional repressor which downregulates the expression of ectodermal genes within the mesoderm. The nuclear protein is structurally similar to the Drosophila snail protein, and is also thought to be critical for mesoderm formation in the developing embryo.

## Application Notes

Titration of the SNAI1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this SNAI1 antibody.

## Storage

Aliquot the SNAI1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.