

## SNAIL Antibody Biotin Conjugate (R35656BTN)

Catalog No.	Formulation	Size
R35656BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Mouse
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6615
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SNAIL antibody is available for research use only.



Western blot testing of A431 lysate with biotinylated SNAIL antibody at 0.1ug/ml.  
Predicted molecular weight ~29kDa.

## Description

Additional name(s) for this target protein: SNAI1, SLUGH2, SNA, SNAH

## Application Notes

Optimal dilution of the SNAIL antibody should be determined by the researcher.

## Immunogen

Amino acids RKPSDPNRKPNY were used as the immunogen for this SNAIL antibody.

## **Storage**

Aliquot and store the SNAIL antibody at -20oC.