


SNTA1 Antibody (R33173)

Catalog No.	Formulation	Size
R33173-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6640
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SNTA1 antibody is available for research use only.



Western blot testing of human muscle lysate with SNTA1 antibody at 0.01ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~54kDa.

Description

Additional name(s) for this target protein: Syntrophin, alpha 1

Application Notes

Optimal dilution of the SNTA1 antibody should be determined by the researcher.

Immunogen

Amino acids ASGRRAPRTGLLE were used as the immunogen for this SNTA1 antibody.

Storage

Aliquot and store the SNTA1 antibody at -20oC.