


SNX12 Antibody (R33140)

Catalog No.	Formulation	Size
R33140-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29934
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SNX12 antibody is available for research use only.



IHC testing of FFPE human breast with SNX12 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Sorting nexin 12

Application Notes

Optimal dilution of the SNX12 antibody should be determined by the researcher.

Immunogen

Amino acids AIDRNYVPGKVRQ were used as the immunogen for this SNX12 antibody.

Storage

Aliquot and store the SNX12 antibody at -20oC.