

## SNX15 Antibody (R33747)

Catalog No.	Formulation	Size
R33747-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29907
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This SNX15 antibody is available for research use only.



Western blot testing of 293 lysate with SNX15 antibody at 0.5ug/ml. Predicted molecular weight: ~38/29kDa (isoforms A/B), observed here at ~45kDa.

## Description

Additional name(s) for this target protein: Sorting nexin 15

## Application Notes

Optimal dilution of the SNX15 antibody should be determined by the researcher.

## Immunogen

Amino acids EEILRLHLSQLPP were used as the immunogen for this SNX15 antibody. This sequence is common to

isoforms A and B.

## **Storage**

Aliquot and store the SNX15 antibody at -20oC.