


SNX18 Antibody (R34685)

Catalog No.	Formulation	Size
R34685-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	112574
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SNX18 antibody is available for research use only.



Western blot testing of A549 lysate with SNX18 antibody at 1ug/ml. Predicted molecular weight: ~72kDa.

Description

Additional name(s) for this target protein: Sorting nexin 18, Sorting nexin associated golgi protein 1, SNAG1

Application Notes

Optimal dilution of the SNX18 antibody should be determined by the researcher.

Immunogen

Amino acids KVQRVPLMTVLSF were used as the immunogen for this SNX18 antibody.

Storage

Aliquot and store the SNX18 antibody at -20oC.