


SNX19 Antibody (R36309)

Catalog No.	Formulation	Size
R36309-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	399979
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This SNX19 antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with SNX19 antibody at 0.3ug/ml.
Predicted/observed molecular weight: ~110kDa.

Description

Additional name(s) for this target protein: Sorting nexin 19

Application Notes

Optimal dilution of the SNX19 antibody should be determined by the researcher.

Immunogen

Amino acids TTSASDTPGNSKR were used as the immunogen for this SNX19 antibody.

Storage

Aliquot and store the SNX19 antibody at -20oC.