


SNX5 Antibody (Biotin Conjugate) (R35966BTN)

Catalog No.	Formulation	Size
R35966BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27131
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SNX5 antibody is available for research use only.



Western blot testing of A549 lysate with biotinylated SNX5 antibody at 2ug/ml.
Predicted/observed molecular weight: ~48kDa.

Description

Additional name(s) for this target protein: Sorting nexin 5

Application Notes

Optimal dilution of the SNX5 antibody should be determined by the researcher.

Immunogen

Amino acids SLLQSCIDLFKNN were used as the immunogen for this SNX5 antibody.

Storage

Aliquot and store the SNX5 antibody at -20oC.