


SNX5 Antibody (R35966)


Catalog No.	Formulation	Size
R35966-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27131
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SNX5 antibody is available for research use only.



Western blot testing of Jurkat lysate with SNX5 antibody at 1ug/ml. Predicted/observed molecular weight: ~48 kDa.



IHC staining of FFPE human breast with SNX5 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Sorting nexin 5

Application Notes

Optimal dilution of the SNX5 antibody should be determined by the researcher.

Immunogen

Amino acids SLLQSCIDLFKNN were used as the immunogen for this SNX5 antibody.

Storage

Aliquot and store the SNX5 antibody at -20oC.