


## SNX8 Antibody (R36242)

Catalog No.	Formulation	Size
R36242-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29886
Applications	Immunohistochemistry (FFPE) : 4-8ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This SNX8 antibody is available for research use only.



IHC testing of FFPE human kidney with SNX8 antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Sorting nexin 8

### Application Notes

Optimal dilution of the SNX8 antibody should be determined by the researcher.

### Immunogen

Amino acids TPPCSPPEDGLCPH were used as the immunogen for this SNX8 antibody.

## Storage

Aliquot and store the SNX8 antibody at -20oC.