


SOCS4 Antibody (R35948)

Catalog No.	Formulation	Size
R35948-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	122809
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SOCS4 antibody is available for research use only.



IHC testing of FFPE human pancreas with SOCS4 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Suppressor of cytokine signaling 4

Application Notes

Optimal dilution of the SOCS4 antibody should be determined by the researcher.

Immunogen

Amino acids SKVRVLRIDAPEQQ were used as the immunogen for this SOCS4 antibody.

Storage

Aliquot and store the SOCS4 antibody at -20oC.