


SOCS7 Antibody (R33377)

Catalog No.	Formulation	Size
R33377-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	30837
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This SOCS7 antibody is available for research use only.



Western blot testing of Daudi lysate with SOCS7 antibody at 1ug/ml. The expected ~60kDa band and the additional ~28kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Suppressor of cytokine signaling 7

Application Notes

Optimal dilution of the SOCS7 antibody should be determined by the researcher.

Immunogen

Amino acids DAGRKPKLTRTQS were used as the immunogen for this SOCS7 antibody.

Storage

Aliquot and store the SOCS7 antibody at -20oC.