


SOX11 Antibody (R36049)

Catalog No.	Formulation	Size
R36049-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6664
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This SOX11 antibody is available for research use only.



Western blot testing of MOLT4 lysate with SOX11 antibody at 0.5ug/ml.
Predicted/observed molecular weight ~47kDa/60kDa (PMID: 19880779).

Description

Additional name(s) for this target protein: SRY (sex determining region Y)-box 11

Application Notes

Optimal dilution of the SOX11 antibody should be determined by the researcher.

Immunogen

Amino acids SRLYYSFKNITKQHP were used as the immunogen for this SOX11 antibody.

Storage

Aliquot and store the SOX11 antibody at -20oC.