


SP140 Antibody (R33784)

Catalog No.	Formulation	Size
R33784-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11262
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SP140 antibody is available for research use only.



Western blot testing of human spleen lysate with SP140 antibody at 2ug/ml. Expected molecular weight: 100/140kDa.

Description

Additional name(s) for this target protein: Speckled 140 kDa, LYSP100-A, LYSP100-B

Application Notes

Optimal dilution of the SP140 antibody should be determined by the researcher.

Immunogen

Amino acids EGQNLQEQVCPEP were used as the immunogen for this SP140 antibody.

Storage

Aliquot and store the SP140 antibody at -20oC.