


SPHK1 Antibody (R33664)

Catalog No.	Formulation	Size
R33664-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8877
Applications	Western Blot : 2-3ug/ml Immunoprecipitation : suitable ELISA (peptide) LOD : 1:16000
Limitations	This SPHK1 antibody is available for research use only.



Western blot testing of human peripheral blood mononucleocyte lysate with SPHK1 antibody at 2ug/ml. Predicted molecular weight: ~43/51/44kDa (isoforms 1/2/3).

Description

Additional name(s) for this target protein: Sphingosine kinase 1

Application Notes

Optimal dilution of the SPHK1 antibody should be determined by the researcher.

Immunogen

Amino acids DVDLESEKYRRLGE were used as the immunogen for this SPHK1 antibody.

Storage

Aliquot and store the SPHK1 antibody at -20oC.