


SPINT2 Antibody (R35896)

Catalog No.	Formulation	Size
R35896-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10653
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This SPINT2 antibody is available for research use only.



Western blot testing of human placenta lysate with SPINT2 antibody at 1ug/ml. Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: Serine peptidase inhibitor, Kunitz type, 2

Application Notes

Optimal dilution of the SPINT2 antibody should be determined by the researcher.

Immunogen

Amino acids SGDDKEQLVKNT were used as the immunogen for this SPINT2 antibody.

Storage

Aliquot and store the SPINT2 antibody at -20oC.