


SRB1 Antibody (R34804)

Catalog No.	Formulation	Size
R34804-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	949
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This SRB1 antibody is available for research use only.



Western blot testing of human adrenal gland lysate with SRB1 antibody at 0.03ug/ml.
Predicted molecular weight: ~57/82kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Scavenger receptor class B, member 1, SR-B1; SCARB1

Application Notes

Optimal dilution of the SRB1 antibody should be determined by the researcher.

Immunogen

Amino acids NGLSKVDFWHSDQ were used as the immunogen for this SRB1 antibody.

Storage

Aliquot and store the SRB1 antibody at -20oC.

References (1)