


## SRD5A1 Antibody (R33183)

Catalog No.	Formulation	Size
R33183-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	6715
<b>Applications</b>	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:64000
<b>Limitations</b>	This SRD5A1 antibody is available for research use only.



IHC testing of FFPE human prostate with SRD5A1 antibody at 4ug/ml. HIER: steamed with pH9 Tris/EDTA buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Steroid-5-alpha-reductase, alpha polypeptide 1

### Application Notes

Optimal dilution of the SRD5A1 antibody should be determined by the researcher.

### Immunogen

Amino acids ATATGVAEERLLC were used as the immunogen for this SRD5A1 antibody.

## **Storage**

Aliquot and store the SRD5A1 antibody at -20oC.