


SRD5A2 Antibody (R34532)

Catalog No.	Formulation	Size
R34532-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6716
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This SRD5A2 antibody is available for research use only.



Western blot testing of human testis lysate with SRD5A2 antibody at 1ug/ml. Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Steroid 5 alpha-reductase 2

Application Notes

Optimal dilution of the SRD5A2 antibody should be determined by the researcher.

Immunogen

Amino acids KMFEDYPKSRK were used as the immunogen for this SRD5A2 antibody.

Storage

Aliquot and store the SRD5A2 antibody at -20oC.