


Stabilin 1 Antibody (R33443)

Catalog No.	Formulation	Size
R33443-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23166
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Stabilin 1 antibody is available for research use only.



IHC testing of FFPE human spleen with Stabilin 1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: STAB1

Application Notes

Optimal dilution of the Stabilin 1 antibody should be determined by the researcher.

Immunogen

Amino acids DFPDTQRILTVK were used as the immunogen for this Stabilin 1 antibody.

Storage

Aliquot and store the Stabilin 1 antibody at -20oC.