


# STAG2 Antibody (R35597)

Catalog No.	Formulation	Size
R35597-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10735
Applications	Western Blot : 0.01ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This STAG2 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: Stromal antigen 2; SA2

### **Application Notes**

Optimal dilution of the STAG2 antibody should be determined by the researcher.

1. This STAG2 antibody will detect isoform a/b.

## Immunogen

Amino acids REPKRLRPEDS were used as the immunogen for this STAG2 antibody.

#### **Storage**

Aliquot and store the STAG2 antibody at -20oC.